**OPERATING INSTRUCTIONS - MODE D'EMPLOI
MANUAL DE INSTRUCCIONES**

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WARNING

N.B.: This "WARNING" is an integral part of the instrument.

This instrument is to be used **STRICTLY** for scientific purposes and/or for internal factory control **ONLY**.

Not for the uses referred to in Art. 1 (2) (a) - L. 90/384 EEC. Except for the models of the CE approved version (legal metrology).

THIS INSTRUMENT IS NOT TO BE USED IN AREAS WHERE A DANGER OF EXPLOSION EXISTS.

Electronic balances, which are sensitive to changes of gravity, must be checked and calibrated in the place where they are to be used.

ALL OUR INSTRUMENTS MEET INTERNATIONAL EMC-EMI-RFI standards according to 89/336 EEC Directive on Electromagnetic Compatibility.

Some models are CE homologated according to 90/384 EEC Directive (legal metrology) and may be used for the manufacture of pre-fabricated packagings, preparation of medical prescriptions at a chemist's, etc.

It is important to note the following points for correct operation of the instrument:

- 1 Read the technical instructions given in the operating manual carefully. This manual contains information about installation and use. If these instructions are not carried out carefully, the accuracy of results may be impaired. If in doubt, please consult the supplier's engineers.
- 2 As for any other electrical equipment, installation should be carried out by a competent person. In particular:
 - ensure that the instrument is correctly earthed, if it is provided for;
 - do not install the instrument in areas of high fire risk, for example, in the presence of inflammable gases and vapours;
 - switch off and disconnect the power supply before removing any cover;
 - do not touch the instrument with wet hands.
- 3 Install the instrument on a vibration free base and away from draughts or sources of heat. Check the instrument with test weightings after installation or after repositioning.
- 4 Avoid connection to a power supply showing high voltage variations due to other loads. The instrument should not be connected to the same electrical circuit as other high power consuming equipment as large fluctuations may cause some inaccuracy.
- 5 The instrument should be checked every 6 months for correct functioning calibration and accuracy of the internal calibration mass by a suitable qualified person.
- 6 A routine calibration check should be carried out each day before use, using a standard mass.
- 7 Check the zero function before each weighing operation.
- 8 If in doubt about any possible malfunction please call a qualified person. Avoid interference by non-qualified people.
- 9 If the instrument is used for weighing food or drink, ensure that other substances are not weighed on the unit which might cause contamination.
- 10 If the instrument is used for weighing chemicals and other similar substances ensure that it is cleaned and any spillage removed to avoid potentially dangerous reactions.
- 11 Always clean the instrument after use to avoid subsequent inaccuracies and/or damage to the equipment.
- 12 Gibertini Elettronica s.r.l. does not accept any liability for damages or other consequences due to tampering of the instrument or its accessories.

ADDITIONAL PROCEDURES FOR WEIGHING SUBSTANCE WITH GREAT ACCURACY

Gibertini Instruments are manufactured to be reliable in use with high accuracy. However, reduced accuracy is always possible due to ambient conditions, incorrect installation or improper use. Where greater accuracy is required, the following procedures should be observed:

- use modern instruments (manufactured within three years) which are fully checked;
- use certified calibration masses before each weighing operation;
- make at least two measurements;
- consider the use of two instruments installed together and make weightings on both units, otherwise use calibrated masses with a certificate of calibration from a recognised institute.

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* Also in the CE metrologic version according to EEC 90/384

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For CE models: the removal of any seals by an unauthorised person will invalidate the product certification.

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WARNING! If the instrument STOPS WORKING because of incorrect operations on the keyboard, DISCONNECT THE MAIN PLUG, wait for some seconds then RECONNECT AGAIN pressing simultaneously ON/OFF key (2) until the display shows "Pnor".

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TECHNICAL FEATURES

- Power source: 220 V -15% / +10% by external power supply, 50 Hz
- Power consumption: 14 VA
- Stabilization time: 2/3 seconds
- 7 digits display, small decimal capacity *(not available for EU600 models)*
- Bar-graph percent used capacity *(not available for EU600 models)*
- ON/OFF display
- Linearity: ± 2 LSD
- Interface: RS232 I/O
- Working temperature: 10-40°C (optimum 15-30°C obligatory for CE version models)
- Balance dimensions: 210x355x120 mm
- Weight: 6,4 kg (7,1 kg EU500)
- Pan diameter:
190 mm (mod. EU8000, EU6000 and EU3000)
150 mm (mod. EU4000 AR, EU1700, EU1000, EU600)
110 mm (mod. EU500)
- Facility for weighing underneath the balance *(not available for EU600 models)*
- Glass protection for EU500.

ACCESSORIES ON REQUEST

- Calibrated mass 1 kg OIML F1 class - 500 g for EU500
- Calibrated masses OIML F1 - E2 class with or without certificate
- Special pans for powders, animals, etc.
- Devices to measure solid and liquid density with printer (PSS and PSL)
- Weighing chamber with 3 sliding doors for EU500.

In addition it is possible to fit:

Controller for variation of display intensity - of particular use in poorly lit environments;

Modification to the integration time for display of the weight - of particular use when weighing animals.

WARNING

Both balance pan and sample to be weighed must be perfectly clean and at the same temperature.

Remember that glass, chinaware and particularly synthetic materials can become charged with static electricity.

This may have an adverse effect on the normal balance performance.

BALANCE DESCRIPTION (Fig. 1)

- 1 BALANCE PAN
- 2 ON/OFF display key
- 3 MODE key
- 4 TARE key
- 5 Display
- 6 Bar-graph
- 7 Hole for suspended weighing under the balance
- 8 Power supply
- 9 RS 232 interface
- 10 Power connection
- 11 Knob and packing bracket
- 12 Level adjusting feet for the models version CE (90/384)
- 13 Spirit level for the models version CE (90/384).

INSTRUCTION TO REMOVE THE LOCKING DEVICE OF THE BALANCE (Fig. 2)

Unscrew the knob on the locking bracket (11) as indicated by the arrow (anti clockwise).

Fit the pan. Retain the locking bracket in case the balance has to be returned to the supplier. This bracket can be refastened to the balance by replacing the knob.

INSTRUCTION FOR WEIGHING UNDER THE BALANCE (Fig. 3)

The balance is supplied with this particular feature to allow use of, for example, the PSS (device for determination of solid density).

Place the instrument on a suitable support and remove the stopper (7) with a screw driver.

It's now possible to see a hook through the hole. Hang the object to be weighed from the hook.

OPERATING INSTRUCTION

INSTALLATION

Place the balance on a rigid, horizontal and vibration-free support, away from air draughts and heat sources. The models in the CE metrologic version have their built-in spirit level.


Connect the instrument (10) to the power supply (8). The flat side of the socket must be upwards.

Connect the power supply to a 220 V line.


The balance is now ON and the display shows a message formed by a letter followed by numbers.

This indicates the instrument program version. This is useful information for our after-sales service.

After about 2 seconds, the display shows:

A rectangular display box containing the text "HELLO" in a monospaced font, followed by a small square icon.

The instrument then makes a display efficiency test as segments start flashing.

A rectangular display box showing a bar graph on the left (four segments lit) and the text "-888888g" in a monospaced font.

After the automatic zerosetting

A rectangular display box showing a bar graph on the left (four segments lit) and the text "0.0" in a monospaced font.

the balance is ready to use.


Before starting operations, allow the instrument a 15 minute warm up time.

In case of urgency, press the TARE key (4) and put the object on the pan. Repeat this operation several times and verify that the value is always the same.

Always leave the instrument ON during the working day. By using the ON/OFF key (2) it's possible to turn OFF the display. In this way the instrument is always ready to use, without the warm up period.

Except numeric values, the display can show two kinds of message: a steady message regarding a normal situation or a flashing message indicating an abnormal situation or requesting intervention by the use.


When an object is put on the pan, the display shows its mass value. If this value is higher than the capacity of the balance, the display shows a flashing message.

A rectangular display box showing a bar graph on the left (four segments lit) and the text "-HI-" in a monospaced font.

WARNING: never put object on the pan with a mass 50% higher than the capacity of the balance.

Also if the mass value is lower than balance capacity, the object should always be put on the pan carefully.

If for any reason the pan is not correctly positioned, the display shows a flashing message.

A rectangular display box containing the text "-LO-" in a monospaced font.

TARE key

With this key (4) it is possible to zeroset all values, positive or negative, until the capacity of the balance is reached.

This key is particularly useful when weighing powders or liquids or an object that requires a container.


Put the empty container on the pan and press the TARE key. Remove it from the balance and fill it.

During this operation, the display shows a negative value equivalent to the container mass.

When the filled container is replaced on the pan, the display will show the net weight.

Europe 4000 AR


When an object is placed on the pan with a mass higher than 2000 g

A rectangular display box showing a bar graph on the left (four segments lit) and the text "2351g" in a monospaced font.

if the TARE key is then pressed

A rectangular display box showing a bar graph on the left (four segments lit) and the text "0.00" in a monospaced font.

the instrument will automatically return to the more accurate range. It is possible to weigh with this accuracy up to the capacity of the balance.

A rectangular display box showing a bar graph on the left (four segments lit) and the text "1585.96g" in a monospaced font.

BAR-GRAPH

The capacity percentage display (Bar-graph) (6) is on the left side of the panel. Each line represents about 10% of the capacity. This information indicates the percent used capacity when the pan is loaded and the balance has just been zeroset.

Calibration

Electronic balances, which are sensitive to changes of gravity, must be checked and calibrated in the place where they are to be used according to EEC Directive 90/384 by means of calibrated masses after installation.

The calibration may be periodically checked using calibrated masses, delivered upon request.

Proceed as the following:

- 1 Press the MODE key (3) until the display shows the "CAL" message and then press the ON/OFF key (2). The display shows:

CAL

until instrument stability is reached.

Then the flashing message appears on the display

LOAD

- 2 Put the calibration mass on the pan. While the instrument is storing the mass value, the display shows the message

CALC

and then the mass value.

The calibration mass value in F1 class (E2 for mod. EU500) must be equal or very near to the maximum capacity of the balance.

If you do not have suitable individual mass, you can temporarily use lower values up to the maximum capacity of the balance.

Ex. 500 g or 100 g, 200 g, etc.

6000 g or 1000 g, 2000 g, etc.

If the calibration mass value is not correct, the display shows the flashing message

Error

The same message is also shown when a mass is not put on the pan within 10 seconds from the "LOAD" message.

In both situation press the TARE key (4). The instrument retains the last calibration mass value.

FOR SPECIAL MODELS EUROPE CAL

Calibration

- 2 Load the calibration mass, turning the knob towards AUTO CAL.

While the instrument is getting the weight value, on the display shows

CALC

followed by the weight value as soon as the instrument has completed the calibration.

Turn the knob back to its initial position.

MONITORING DATA

(only for EUROPE 500 CAL)

Connecting the computer, it is possible to get on the monitor 4 decimal figures displayed

0.0000

ADDITIONAL FUNCTION: PIECE-COUNTING (CPZ)

(not available for EU600 models)

If you need to use a container, place it on the pan and press TARE key (4).

Press MODE key (3) until the display shows

CPZ

and then ON/OFF key (2) to confirm.

The message

P 10

appears on the display to indicate the number of pieces of the sample.

It is possible to select samples with 10, 20, 50 or 100 pieces.

If the pieces are not uniform, it is recommended that a large number of samples are checked.

Press MODE key (3) until the display shows the desired sample number.

Add the required pieces of the sample and then press ON/OFF key (2) to confirm.

The message is now flashing. Press ON/OFF key again to store the sample weight. It is now possible to add other pieces up to the instrument capacity. The display shows the number of pieces on the pan.

If the weight of the pieces is required, press MODE key (3) until the display shows:

W 10

and then press ON/OFF key (2) to confirm. The display shows the weight value. To return to CPZ function press MODE key until the display shows the corresponding message and then press ON/OFF key to confirm and then TARE key.

Warning! The individual sample weight cannot be lower than the balance accuracy.

If it is lower, the display shows the blinking message:

Error

Press ON/OFF key to cancel. The display shows the last stored sample.

MESSAGE MENU

Steady messages don't require any intervention.

Flashing messages indicate an irregular situation or request of intervention by the operator.

Initial message

Revision level of the internal program (Firmware)

---- Zerosetting operating. Press "T" key.

Er18 Unstable value during zerosetting operation. Press "T" to cancel. The instrument retains the last zero. If the message appears when switching ON, you cannot cancel. Switch OFF for a few seconds and then ON again. If the display shows the same message, the instrument is out of order.

Er19 Unstable value during the calibration. Press "T" to cancel. The instrument keeps the last calibration mass value.

Er23 Loaded pan at switching ON. Press "T" to cancel. The weight placed on the pan is considered as an initial tare.

ErrC Calibration error. The value of the calibration mass is not correct or operating error. Press "T" to cancel. The instrument retains the last calibration mass value.

- **Hi** - The mass is higher than the balance capacity. Remove the objects.

- **Lo** - The mass is lower than the acceptable minimum. The pan is not present or not correctly placed.

ACAL Calibration procedure starts. The message is shown until stability is reached.

LOAD The instrument is waiting for calibration mass.

CALC The instrument is calculating the calibration factor. At the end the display shows the mass value or the "ErrC" message.

RS232 OUTPUT

Host computer communications

Communications between the Balance and Computer can be bidirectional or from Balance to Computer only. The communications parameters can be selected from the following:

- 1 - Baud rate 75 - 9600 bit/sec.
- 2 - Parity Even, Odd or None.
- 3 - Word length 7 or 8 bits.
- 4 - Continuous or On request modes.

The default parameters are 9600 bit/sec, Even Parity, 7 bits/word, 1 stop bit, On request mode.

If different parameters are requested, please, contact our technical service.

RS232 SERIAL INTERFACE

The instrument has a RS232 serial interface (9), sub DB9 type female connector.

9 Pin connection table

Pin	Description	Direction
1	Common of the serial data send command switch	-
2	Transmit Data	0
3	Receive Data	1
4	Data Set Ready	1
5	Signal Ground	-
6	Data Terminal Ready	0
7	N.C.	-
8	Request to Send	-
9	Serial data send command switch	-

The DTR (pin 6) and RTS (pin 8) signals are always active.

The DSR (pin 4) signal is not active.

The serial data send command is a N.O. switch connected between pins 9 and 1.

CONNECTION CABLES TO HOST COMPUTER

If the connector is a 9 pin one, the cable should be as follows:

Computer (9 pin female connector)	EUROPE (9 pin male connector)
2	2
3	3
5	5

If the connector is a 25 pin one, the cable should be as follows:

Computer (25 pin female connector)	EUROPE (9 pin male connector)
3	2
2	3
7	5

Warning! If the instrument does not work because of incorrect operations on the keyboard, disconnect the main plug, wait for some seconds and then reconnect pressing simultaneously ON/OFF key (2) until the display shows "Pnor".

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
Fig. 2

UK Device for determination of the density of solids with balances EUROPE series

- Put the balance on the plexiglass support and remove the protective stopper of the hook. If necessary use a screwdriver (Fig.1).
- Hang the basket "B" to the hook and place the elevator and the glass container (Fig. 2).
- Fill the container with distilled water at 20°C and adjust the elevator so that the water level cover the basket wire knot 2÷5 mm at least. Ensuring that the basket doesn't touch the container. Then zero the balance.
- Weight the solid sample in air (on the pan "A"), note the weight and then ZERO THE BALANCE.
- Transfer the sample into the water (in the basket "B") and read the negative weight (thrust) shown on the display. Ensuring that there are not air bubbles on the basket and on the sample surface.
- The ratio between the weight in air and the thrust is the density of the sample.

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
UK Device for determination of the density of liquids with balances Mod. EUROPE 500

- Remove the the glass cylinder cover and, carefully, the weighing pan.
- Mount the PSL device, like the picture, ensuring that the float support doesn't touch the cylinder support.
- Suspend the float (20 ml at 20°C), and then ZERO THE BALANCE.
- Fill the cylinder with 80 ml of liquid at 20°C; introduce the float and carefully shake it to, reach the liquid temperature. Pay attention that there are not air bubbles on the float surface.
- Suspend again the float on the support and note the negative value shown on the display. For example 19,387g. Pay attention that the liquid cover the basket wire knot 2÷5 mm at least.
- To obtain the liquid density divide the found value for the float volume (20 ml).
- **THE VOLUME CALIBRATION OF THE FLOATER MUST BE PERIODICALLY CHECKED USING DISTILLED WATER AT 20°C. NEVER TOUCH THE GLASS FLOATER WITH THE FINGERS AND KEEP IT ALWAYS VERY CLEAN.**

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